



# Expression Systems Bulk Packaging

## Bulk Liquid Packages

Our "Media Delivery Bags" have the easy-to-use connectors that makes them a snap to connect with your system.

Fill volumes range in size from 4 liters to 200 liters or greater.

Standard tubing and fittings configuration for 50 Liter bags:

082 clear C-Flex 1/4 x 1/8" with Luer Locks

082 clear C-Flex 3/8 x 5/8" with Male Quick Connect


082 clear C-Flex 1/4 x 7/16" with plug

## 50 Liter Drum



## Easy to use bags

10 and 20 liter bags are equipped with hanging bar, one male and one female Luer Locks, male Quick Connect and Injection Port



A - Injection Port

B - Quick Connect

C - Male and Female Luer Lock

Tubing available in clear 082 C-Flex 1/8 X 1/4" with male and female Luer Lock or the larger 3/8 X 5/8" with male Quick Connect.



EXPRESSION SYSTEMS

26 Harter Ave, Suite 2 • Woodland, CA • 95776

Phone: 530.662.2946 • fax 530.662.5359

www.expressionsystems.com