



Expression  
Systems

Doc ID: 15802

Rev. No: 4

## Certificate of Analysis

**Product:** ESF 921 Insect Cell Culture Medium, Protein Free

**Catalog Number:** 96-001

**Lot Number:** 21708039

**Manufacture Date:** 21AUG17

Test	Specification	Result
Appearance	Clear, yellow solution	Pass
Sterility – current USP	Negative for microbial and fungal growth	Pass
pH	$6.10 \leq \text{pH} \leq 6.25$	6.18
Osmolality	400-500 mOsm/kg H <sub>2</sub> O	443 mOsm/kg
Growth on Sf9 and Tni	Pass	Pass

Product is intended for Research Use Only. Not for human use.

This material has been evaluated by Quality Assurance and is approved.

Date: 11SEP17

Donald Le  
Quality Assurance Specialist  
Expression Systems LLC

Expression Systems LLC  
2537 2nd Street, Davis, California 95618  
P: (530)747-2035 F: (530)747-2034  
expressionsystems.com