



An **Advancion** Company

Certificate of Analysis

Product: ESF 921 Insect Cell Culture Medium, Protein Free

Catalog Number: 96-001

Lot Number: 22502014

Manufacture Date: 03MAR25

Test	Specification	Result
Appearance	Clear, yellow solution	Pass
Sterility – USP <71>	Negative for microbial and fungal growth	Pass
pH – USP <791> @25C	6.00 ≤ pH ≤ 6.40	6.17 pH
Osmolality – USP <785>	400-500 mOsmol/kg	453 mOsmol/kg
Endotoxin	< 20 EU/mL	≤0.453 EU/mL
Growth on Sf9 and Tni	Pass	Pass

For Research Use or Further Manufacturing.
Not for diagnostic use or direct administration into humans or animals.

This material has been evaluated by Quality Assurance and is approved.

A handwritten signature in black ink, appearing to read "Lorena Soto".

Date: 18 Mar 25

Lorena Soto De Almanza
Quality Assurance Specialist III
Expression Systems LLC

Expression Systems LLC
2537 2nd Street, Davis, California 95618
P: (530)747-2035 F: (530)747-2034
expressionsystems.com