



An **Advancion** Company

Certificate of Analysis

Product: ESF AF Insect Cell Culture Medium, Animal Free

Catalog Number: 99-300


Lot Number: 22505051

Manufacture Date: 20MAY25

Test	Specification	Result
Appearance	Clear, light brown solution	Pass
Sterility – USP<71>	Negative for microbial and fungal growth	Pass
pH – USP <791> @25C	6.00 ≤ pH ≤ 6.40	6.24 pH
Osmolality – USP <785>	400-500 mOsmol/kg	465 mOsmol/kg
Endotoxin	< 20 EU/mL	<0.400 EU/mL
Growth on Sf9 and Trn	Pass	Pass

For Research Use or Further Manufacturing.
Not for diagnostic use or direct administration into humans or animals.

This material has been evaluated by Quality Assurance and is approved.



Date: 10 Jun 25

Lorena Soto
Quality Assurance Specialist III
Expression Systems LLC

Expression Systems LLC
2537 2nd Street, Davis, California 95618
P: (530)747-2035 F: (530)747-2034
expressionsystems.com