



An **Advancion** Company

Certificate of Analysis

Product: Tni Cells, Frozen in ESF 921, 50M cells/vial

Catalog Number: 94-002F

Lot Number: 22602016

Manufacture Date: 19Feb26

Organism: *Trichoplusia ni*
Identity: Insect
Strain: Tni
Volume: 1 mL/vial
Culture at Time of Harvest: Density- 7.353×10^6 vc/mL
Viability-98.0 %,
Avg cell diameter-19.93 μ m
Storage Medium: 90% ESF 921, 10% DMSO, 0.15M Trehalose

Test	Specification	Result
Gross Sterility (Visual)	Negative for Bacterial and Fungal Growth	Pass
Viability Upon Thaw	> 90%	Pass
Mycoplasma (PCR)	Negative for Mycoplasma	Pass

This material has been evaluated by Quality Assurance and is approved.

A handwritten signature in black ink, appearing to read "Carina Williams", written over a horizontal line.

Carina Williams
Quality Assurance Specialist III
Expression Systems LLC

Date: 14 APR 26

Expression Systems LLC
2537 2nd Street, Davis, California 95618
P: (530)747-2035 F: (530)747-2034
expressionsystems.com